

# LabSpion: Lamp Bracket Model II

The new LabSpion Lamp Bracket Model II is available as an accessory for all existing LabSpion systems. The Model II bracket makes it even easier to attach all sorts of light sources and luminaires to the goniometer.

## Advantages

- The detailed hole pattern creates many possibilities for mounting light sources
- Fits all LabSpion models
- Detach the whole bracket and take it to the workshop
- Fast and easy attachment to the gonio tower with two sturdy hooks and two large thumb screws




Replace your existing lamp bracket with the new base plate



Hook on your new lamp bracket and tighten the two thumb screws



Having really big light sources to measure? Get the extender set.




Fits all LabSpion versions



# SPECIFICATIONS

Detach and remove




**Lamp Bracket Model II**  
**Physical specifications**

Contents:  
Lamp bracket incl. c-plane mount, angle brackets (2 pcs.), fastener handles M8 (2 pcs.)

Dimensions (L x W x H)  
Weight  
Middle hole diameter  
Long slots, width  
Smaller slots, width  
Materials

810 x 150 x 8 mm  
Bracket 2.02 kg + 1.36 kg (c-plane mount) + 0.55 kg (angle brackets and handles)  
70.0 mm  
10.2 mm  
8.2 mm  
Bracket: Powder coated aluminium  
C-plane mount: Powder coated steel and aluminium



**Lamp bracket extension set**  
**Physical specifications**

Contents:  
Lamp bracket extensions (2 pcs.), thumb screws (8 pcs.)

Dimensions (L x W x H)  
Weight  
Middle hole diameter  
Long slots, width  
Smaller slots, width  
Materials

720 x 70 x 8 mm  
0.54 kg + 0.15 kg (thumb screws)  
70.0 mm  
10.2 mm  
8.2 mm  
Brackets: Powder coated aluminium

